SUPPLÉMENT DE LIQUIDE SYNOVIAL

CONTENU

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SYNOVIAL FLUID SUPPLEMENT

CONTENTS

Active ingredient: Sodium Hyaluronate 25 mg.
Excipients: Disodium Phosphate Dodecahydrate, Sodium Dihydrogen Phosphate Dihydrate, Sodium Chloride, Sterile Water for Injection q.s. to 2.5 ml.
pH 6.8~7.8

DESCRIPTION

HYA**FLEX** SYNOVIAL FLUID SUPPLEMENT, is a product supplied in a glass syringe with hyaluronic acid. HYA**FLEX** is a clear, lucid, viscous, sterile and nonpyrogenic solution. The hyaluronic acid (molecular weight 650000~1200000 Da) is extracted from microorganism. The product is for single use only. Hyaluronic acid is a natural polysaccharide which occurs as an important structural element of skin, subcutaneous and connective tissues, such as synovial fluid. Hyaluronic acid is identical in all living organisms and highly biocompatible. Each one ml of HYA**FLEX** contains 10 mg of sodium hyaluronate dissolved in physiology saline.

PACKAGE

2.5 ml per syringe

INDICATION

HYA**FLEX** is indicated for treatment of pain in osteoarthritis (OA) joint.

USAGE

HYA**FLEX** is used 3 syringes per treatment, 1 syringe per week for adult into the osteoarthritis (OA) joint. HYA**FLEX** is used for injecting to joint. The dosage and injection time is according to symptom of patient. Anticipated duration of efficacy is at least 6 months after HYA**FLEX** treatment. HYA**FLEX** is used with a 21-25 G needle by injecting the material into articular cavity of joint.

THERAPEUTIC EFFICACY

HYA**FLEX** can counterbalance the mechanical stress transmitted from outside. It is a barrier protecting the joint. HYA**FLEX** can protect synovial tissues and promote healing of articular cartilage.

CONTRAINDICATIONS

- Do not administer to patients with known hypersensitivity to hyaluronate preparations.
- Intra-articular injections are contraindicated in cases of past and present infections or skin diseases in the area of the injection site.

WARNING

HYA**FLEX** is intended for use as intra-articular implant. Don't re-sterilize HYA**FLEX** or mix with other products. Injections into blood vessels and surrounding tissue should be avoided. Doctor should take care of patients who ever got allergy about other drugs and who have liver disease. If re-using, it may cause infection even seriously inflammatory response.

USAGE FOR OLD MEN

Because the organ functions of old men are weak, doctors should pay attention to them.

USAGE FOR MOTHER TO BE AND LACTATION

If anticipated efficacy of treatments are better than adverse events, HYA**FLEX** could be intended to use. During the treatment period, lactation mothers should stop to suckle.

USAGE FOR CHILDREN

The safety of use for children is not sure; therefore, the doctor should pay attention when treating.

ANTICIPATED SIDE-EFFECTS

The major side-effects are pain, swelling and warmth localized to the treated knee.

NOTE!

HYA**FLEX** is only intended to be administered by authorized personnel in accordance with local legislation. Do not over screw the plug rod to prevent product leakage, and avoid stopper deformation that increases the extrusion force. The syringe, the needle and any unused material must be discarded directly after the treatment session. If the package is damaged, do not use. Before use, please refer to instruction.

STORAGE

As indication on package, store under room temperature. Protect from sunlight and freezing.

MANUFACTURE DATE AND EXPIRY DATE

As indicated on product label and package.

SHELF LIFE



3 years

 SciVision Biotech Inc.
9, South 6th Rd., K.E.P.Z., Kaohsiung 806, Taiwan, R.O.C.

EXCLUSIVE DISTRIBUTOR FOR THE BENELUX COUNTRIES

S.A. Laboratoires Pharmaceutiques Trenker
Avenue Thomas Edison 32 • BE 1402 Thines • Belgium

SYMBOLS ON THE PACKAGING

	Do not re-use		Do not use if package is damaged
	Temperature limitation		Lot number
	Date of manufacture		Manufacturer
	Use-by date		Consult instructions for use
	Sterile. The contents of the syringe have been sterilized by using moist heat.		
	CE marked, 2797 is No. of Notified Body.		

 **Obelis s.a,**
Bd. General Wahis 53, 1030 Brussels • Belgium
T +32 (0)2.732.59.54
F +32 (0)2.732.60.03
E mail@obelis.net