Quies® Docuspray

2272/24F
REV1-01/06/2024

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MD Dispositivo medico.

Ⓜ EAR HYGIENE

Regular care of the ear
For easy removal of earwax
and dead skin cells

With marine origin trace elements
Gas propellant-free ear spray
Easy to use, with single dose pump

INDICATIONS

Quies Docuspray cleans the ear gently and prevents earwax blockages from forming thanks to its complete formula, combining:

- Marine origin trace elements
- Sodium bicarbonate
- Polysorbate-80 which promotes earwax release.

RECOMMENDATIONS FOR USE

- Easy to use, with single dose pump.
- Rinse-free. Gas propellant-free spray.
- Ensure that the cone is fixed to the nozzle. Place it at the entrance to the ear canal and press the pump. Massage the base of the ear. Tilt the head to encourage liquid and residues to flow.
- Wipe with a clean tissue. Rinse the nozzle in hot water after each use. Leave it in hot water for 30 minutes if it is blocked.
- 1 to 2 sprays, 2 to 3 times per week.
- 1 to 2 sprays per day, in case of heavy earwax production or a dirty environment.
- Can be used for children aged over 30 months and adults. Suitable for hearing aid wearers, people producing high levels of earwax, after having removed their hearing aids.

COMPOSITION

Sodium bicarbonate: 1,1%, Sea salt: 0,1%, Polysorbate-80: 0,1%, Purified water, Benzalkonium chloride, Mono and disodium phosphate

PRECAUTIONS FOR USE

Keep out of the sight and reach of children. In case of contact with mucosa, rinse with clean water. If in doubt about using the device, ask your doctor for advice. If you experience a serious incident related to this device, contact quies@quies.com and inform the appropriate authority.

CONTRAINDICATIONS

Do not use if you have an ear infection, perforated eardrum or are wearing tympanostomy tubes or after ear surgery. Do not use if you are known to be allergic to one of the ingredients. Do not use on damaged skin.

STORAGE:

It is recommended to use the product within the year of opening the bottle. Any used or started device must be destroyed following local regulations regarding waste disposal.

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