Quies® Doculyse

22298/24H
REV.2-01/11/2024

(FR) ANTI-BOUCHON DE CÉRUMEN

Spray auriculaire

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(GB) EARWAX BLOCKAGE REMOVER

Ear spray

Helps dissolve earwax blockages and remove them from the external ear canal. ENT specialist tested. Earwax blockages can alter hearing, cause buzzing or dizziness.

Quies Doculyse allows the earwax blockage to be dissolved and removed thanks to its main active ingredient, Sodium Docusate. Its high tolerance formula has been subject to tests under ENT specialist control.

RECOMMENDATION FOR USE

- Easy to use, with single dose pump.
- Gas propellent-free spray.
- Ensure that the cone is fixed to the nozzle. Place it at the entrance to the ear canal and press the pump. Massage the base of the ear. Leave to act for 15 minutes by tilting the head to the opposite side of the ear being treated. Then, in order to encourage liquid and residues to flow, tilt the head to the other side. Rinse-free.
- Wipe with a clean tissue. Rinse the nozzle in hot water after each use. Leave it in hot water for 30 minutes if it is blocked.
- To remove a recent blockage:
 1 to 2 sprays, repeat it necessary.
- To remove an old blockage:
 1 to 2 sprays, 3 times per day the 3 or 4 days preceding removal by an ENT specialist.
- As prevention: 1 to 2 sprays per week.
- Can be used for babies aged over 6 months and adults. Suitable for people with an earwax blockage or heavy earwax secretion. Suitable for hearing aid wearers, people producing high levels of earwax, after having removed their hearing aid.

COMPOSITION

Sodium Docusate: 0.5%, Sodium citrate, Citric acid, Purified water, Phenoxyethanol, Methylparaben, Ethylparaben, Butylparaben, Propylparaben.

PRECAUTIONS FOR USE

Keep out of the sight and reach of children. In case of pain in the ear, consult your doctor. In case of contact with mucosa, rinse with clean water. If in doubt about using the device, ask your doctor for advice. If you experience a serious incident related to this device, contact quies@quies.com and inform the appropriate authority.

CONTRAINDICATIONS

Do not use if you have an ear infection, perforated eardrum or are wearing tympanostomy tubes or after ear surgery. Do not use in case of irritation. Do not use if you are known to be allergic to one of the ingredients.

STORAGE

It is recommended to use the product within the year of opening the bottle. Any used or started device must be destroyed following local regulations regarding waste disposal.

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